

Product Datasheet

Rat MSH Alpha ELISA Kit (orb410700)

Description	Rat Pro-opiomelanocortin(POMC) ELISA KIT
Reactivity	Rat
Range	7.8 pg/mL-500 pg/mL
Target	POMC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P01194
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.