

## Product Datasheet

### Rat IL17A ELISA Kit (orb410712)

<b>Description</b>	Rat Interleukin 17,IL-17 ELISA KIT
<b>Reactivity</b>	Rat
<b>Range</b>	93.75 pg/mL-6000 pg/mL
<b>Target</b>	IL17A
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>Q61453</b>
<b>Sensitivity</b>	23.43 pg/mL
<b>Expiration Date</b>	Please enquire.