

Product Datasheet

Rat Estradiol ELISA Kit (orb410777)

| | |
|--------------------------|------------------------------|
| Description | Rat Estradiol(E2) ELISA Kit |
| Reactivity | Rat |
| Range | 1.25 ng/mL-50 ng/mL |
| Target | E2 Ab |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma |
| Sensitivity | 1.25 ng/mL |
| Expiration Date | Please enquire. |