

## Product Datasheet

### Rat Oxytocin ELISA Kit (orb410791)

<b>Description</b>	Rat Oxytocin ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	7.5 pg/mL-600 pg/mL
<b>Target</b>	OT
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma.
<b>Sensitivity</b>	9.375 pg/mL
<b>Expiration Date</b>	Please enquire.