

## Product Datasheet

### Rat MCP1 ELISA Kit (orb410795)

<b>Description</b>	Rat monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP-1/MCAF ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	CCL2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P14844</b>
<b>Sensitivity</b>	7.8 pg/mL
<b>Expiration Date</b>	Please enquire.