

## Product Datasheet

### Rat IGFBP5 ELISA Kit (orb410807)

**Description** Rat Insulin-like growth factor binding protein 5,IGFBP-5 ELISA kit

**Reactivity** Rat

**Range** 0.312 ng/mL-20 ng/mL

**Target** IGFBP5

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma,tissue homogenates

**Uniprot ID** [P24594](#)

**Sensitivity** 0.078 ng/mL

**Expiration Date** Please enquire.