

## Product Datasheet

### Rat MMP3 ELISA Kit (orb410935)

**Description** Rat matrix metalloproteinase 3, MMP-3  
ELISA kit

**Reactivity** Rat

**Range** 1.56 pg/mL-100 pg/mL

**Target** MMP3

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell culture supernates,  
tissue homogenates

**Uniprot ID** **P03957**

**Sensitivity** 0.39 pg/mL

**Expiration Date** Please enquire.