

## Product Datasheet

### Rat FGF2 ELISA Kit (orb410952)

<b>Description</b>	Rat basic fibroblast growth factor,bFGF ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312 ng/ml -20 ng/ml
<b>Target</b>	FGF2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P13109</b>
<b>Sensitivity</b>	0.078 ng/ml
<b>Expiration Date</b>	Please enquire.