

Product Datasheet

Rat CCL20 ELISA Kit (orb410976)

| | |
|--------------------------|---|
| Description | Rat macrophage inflammatory protein 3 α ,MIP-3 α ELISA Kit |
| Reactivity | Rat |
| Range | 0.156 ng/ml - 10 ng/ml |
| Target | CCL20 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P97884 |
| Sensitivity | 0.039 ng/ml |
| Expiration Date | Please enquire. |