

Product Datasheet

Rat LOX-1 ELISA Kit (orb411004)

Description Rat Soluble Lectin-like Oxidized Low Density Lipoprotein Receptor-1(sLOX-1)ELISA Kit

Reactivity Rat

Range 0.156 ng/mL-10 ng/mL

Target OLR1

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **070156**

Sensitivity 0.039 ng/mL

Expiration Date Please enquire.