

## Product Datasheet

### Rat BNP ELISA Kit (orb411047)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Rat brain natriuretic peptide, BNP ELISA Kit |
| <b>Reactivity</b>        | Rat  |
| <b>Range</b>             | 125 pg/mL-8000 pg/mL                         |
| <b>Target</b>            | BNP  |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                 |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates            |
| <b>Uniprot ID</b>        | <b>P13205</b>                                |
| <b>Sensitivity</b>       | 31.25 pg/mL                                  |
| <b>Expiration Date</b>   | Please enquire.                              |