

## Product Datasheet

### Rat NTS ELISA Kit (orb411122)

<b>Description</b>	Rat neurotensin (NTS) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Target</b>	NTS
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P20068</b>
<b>Sensitivity</b>	3.9 pg/mL
<b>Expiration Date</b>	Please enquire.