

## Product Datasheet

### Rat SST ELISA Kit (orb411151)

<b>Description</b>	Rat Somatostatin,SS ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	1.56 pg/mL-100 pg/mL
<b>Target</b>	SST
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P60042</b>
<b>Sensitivity</b>	0.39 pg/mL
<b>Expiration Date</b>	Please enquire.