

Product Datasheet

Rat CCL11 ELISA Kit (orb411155)

Description	Rat Eotaxin 1 (Eotaxin 1/CCL11/ECF)ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml -1000 pg/ml
Target	CCL11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P97545
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.