

## Product Datasheet

### Rat LEP ELISA Kit (orb411158)

<b>Description</b>	Rat Leptin,LEP ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.78 ng/ml-50 ng/ml.
<b>Target</b>	LEP
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates.
<b>Uniprot ID</b>	<b>P50596</b>
<b>Sensitivity</b>	0.195 ng/ml.
<b>Expiration Date</b>	Please enquire.