

## Product Datasheet

### Rat RLN1 ELISA Kit (orb411167)

<b>Description</b>	Rat relaxin 1 (RLN1) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	2.4 pg/mL-150 pg/mL
<b>Target</b>	RLN1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P01347</b>
<b>Sensitivity</b>	0.6 pg/mL
<b>Expiration Date</b>	Please enquire.