

Product Datasheet

Rat AGE ELISA Kit (orb411205)

| | |
|--------------------------|--|
| Description | Rat advanced glycation end products,AGEs ELISA Kit |
| Reactivity | Rat |
| Range | 3.9 µg/mL-250 µg/mL |
| Target | AGEs |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Sensitivity | 0.97 µg/mL |
| Expiration Date | Please enquire. |