

Product Datasheet

Rat SHBG ELISA Kit (orb411247)

Description	Rat sex hormone-binding globulin (SHBG) ELISA kit
Reactivity	Rat
Range	375 ng/mL-6000 ng/mL
Target	SHBG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P08689
Sensitivity	375 ng/mL
Expiration Date	Please enquire.