

Product Datasheet

Rat Vasoactive Intestinal Peptide ELISA Kit (orb411252)

Description	Rat Vasoactive Intestinal Peptide, VIP ELISA Kit
Reactivity	Rat
Range	1.56 pg/mL-100 pg/mL
Target	VIP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P32649, P01283
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.