

## Product Datasheet

### Rat Ang1-7 ELISA Kit (orb411255)

<b>Description</b>	Rat angiotensin 1-7 (Ang1-7) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	18.75 pg/mL-1200 pg/mL
<b>Target</b>	Ang1-7
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	4.68 pg/mL
<b>Expiration Date</b>	Please enquire.