

## Product Datasheet

### Rat VWF ELISA Kit (orb411275)

**Description** Rat von Willebrand Factor,vWF ELISA Kit

**Reactivity** Rat

**Range** 0.312 ng/mL-20 ng/mL

**Target** VWF

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [Q62935](#)

**Sensitivity** 0.078 ng/mL

**Expiration Date** Please enquire.