

## Product Datasheet

### Rat PDGF-BB ELISA Kit (orb411295)

<b>Description</b>	Rat platelet-derived growth factor-BB,PDGF-BB ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	39.063 pg/mL-2500 pg/mL
<b>Target</b>	PDGF-BB
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Sensitivity</b>	9.77 pg/mL
<b>Expiration Date</b>	Please enquire.