

Product Datasheet

Rat INSL3 ELISA Kit (orb411311)

Description	Rat Insulin-Like 3 (INSL3) ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	INSL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9WUK0
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.