

Product Datasheet

Rat RANTES/CCL5 ELISA Kit (orb411356)

Description	Rat regulated on activation in normal T-cell expressed and secreted/C-C motif chemokine 5 (RANTES/CC...
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	RANTES/CCL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P50231
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.