

Product Datasheet

Rat Neutrophil Elastase ELISA Kit (orb411365)

Description	Rat neutrophil elastase,NE ELISA Kit
Reactivity	Rat
Range	12.5 ng/mL-800 ng/mL
Target	NE
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P08246
Sensitivity	3.12 ng/mL
Expiration Date	Please enquire.