

## Product Datasheet

### Rat CTSB ELISA Kit (orb411388)

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Description</b>       | Rat Cathepsin B(CTSB) ELISA kit                    |
| <b>Reactivity</b>        | Rat                                                |
| <b>Range</b>             | 18.75 ng/mL-1200 ng/mL                             |
| <b>Target</b>            | CTSB                                               |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell lysates |
| <b>Uniprot ID</b>        | <b>P00787</b>                                      |
| <b>Sensitivity</b>       | 4.67 ng/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                    |