

## Product Datasheet

### Sheep IL8 ELISA Kit (orb411451)

**Description** Sheep Interleukin-8,IL-8 ELISA Kit

**Reactivity** Sheep

**Range** 6.25 pg/mL-400 pg/mL

**Target** IL8

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue  
homogenates

**Uniprot ID** **P36925**

**Sensitivity** 1.56 pg/mL

**Expiration Date** Please enquire.