



Product Datasheet

hnRNP A2B1/HNRNPA2B1 Antibody (orb412956)


Description	hnRNP A2B1/HNRNPA2B1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human hnRNP A2B1 (KTLETVPLERKKREKEQFRKLFIGGLSFETTEESLRNYEQWGKL).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P22626
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A431 cells us...



WB analysis of hnRNP A2B1 using anti-hnR...



IF analysis of hnRNP A2B1 using anti-hnR...