



Product Datasheet

Fatty Acid Binding Protein 5/FABP5 Antibody (orb412979)


Description	Fatty Acid Binding Protein 5/FABP5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide corresponding to a sequence of mouse FABP5 (KWRLMESHGFEEYMKELGVGLALRKMAAMAKPD).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	Q05816
Expiration Date	12 months from date of receipt.



WB analysis of FABP5 using anti-FABP5 an...



IF analysis of FABP5 using anti-FABP5 an...



IHC analysis of FABP5 using anti-FABP5 a...