



Product Datasheet

PLCG 2/PLCG2 Antibody (orb412987)


Description	PLCG 2/PLCG2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E. coli-derived human PLCG 2 recombinant protein (Position: R1201-V1258).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	148 kDa
Uniprot ID	P16885
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of THP-1 cells using anti-PLCG2 antibody.



Western blot analysis of PLCG 2/PLCG2 using anti-PLCG2 antibody.



Immunofluorescence analysis of THP-1 cells showing PLCG 2/PLCG2 expression using anti-PLCG2 antibody.