



Product Datasheet

BCMA/TNFRSF17 Antibody (orb413027)


Description	BCMA/TNFRSF17 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E. coli-derived human BCMA recombinant protein (Position: M1-A54).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	Q02223
Expiration Date	12 months from date of receipt.



WB analysis of BCMA using anti-BCMA anti...



Sandwich ELISA - Recombinant human BCMA/...



IHC analysis of BCMA using anti-BCMA ant...