

## Product Datasheet

### BAR 502 (orb413247)

<b>Description</b>	BAR502 is a dual FXR and GPBAR1 agonist with IC50 values of 2 µM and 0.4 µM, respectively.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>25</sub> H <sub>44</sub> O <sub>3</sub>
<b>MW</b>	392.62
<b>CAS Number</b>	[1612191-86-2]
<b>SMILES</b>	CCC1[C@@H]2CC(CC[C@@]2([C@@H]3CC[C@]4([C@H])([C@@H]3C1O)CC[C@@H]4C(C)CCO)C)C)O
<b>Expiration Date</b>	12 months from date of receipt.