

Product Datasheet

GST antibody (orb413654)

Description	Mouse monoclonal antibody against GST conjugated to PerCP
Species/Host	Mouse
Reactivity	All
Conjugation	PerCP
Tested Applications	ELISA, WB
Immunogen	Glutathione-S-Transferase (GST) tagged protein
Target	GST
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-549 was sufficient for detection of GST in 1 µg of Purified GST tagged PmtA by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	P08515
Dilution Range	WB (1:1000); ELISA (1:800)
Expiration Date	12 months from date of receipt.