
Product Datasheet

**DYKDDDDK Tag antibody (Equivalent to Sigma's
Anti-FLAG M2 antibody) (orb413729)**

Description

Rabbit polyclonal antibody to

Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to the DYKDDDDK Tag, conjugated to Keyhole Limpet Haemocyanin (KLH).
Target	DYKDDDDK Tag
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-720 was sufficient for detection of DYKDDDDK Tag in 15 µg of DYKDDDDK (FLAG) Positive Control Lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.