
Product Datasheet

**DYKDDDDK Tag antibody (Equivalent to Sigma's
Anti-FLAG M2 antibody) (orb413740)**

Description

Rabbit polyclonal antibody

conjugated to Biotin

Species/Host

Rabbit

Reactivity

All

Conjugation

Biotin

Tested Applications

WB

Immunogen

Synthetic peptide corresponding to the DYKDDDDK Tag, conjugated to Keyhole Limpet Haemocyanin (KLH).

Target

DYKDDDDK Tag

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-720 was sufficient for detection of DYKDDDDK Tag in 15 µg of DYKDDDDK (FLAG) Positive Control Lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.