

---

## Product Datasheet

**DYKDDDDK Tag antibody (Equivalent to Sigma's  
Anti-FLAG M2 antibody) (orb413745)**

**Description**

Rabbit polyclonal antibody

conjugated to RPE

**Species/Host**

Rabbit

**Reactivity**

All

**Conjugation**

RPE

**Tested Applications**

WB

**Immunogen**

Synthetic peptide corresponding to the DYKDDDDK Tag, conjugated to Keyhole Limpet Haemocyanin (KLH).

**Target**

DYKDDDDK Tag

**Preservatives**

PBS pH 7.4, 50% glycerol, 0.09% sodium azide \*Storage buffer may change when conjugated

**Concentration**

1 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

A 1:1000 dilution of SPC-720 was sufficient for detection of DYKDDDDK Tag in 15 µg of DYKDDDDK (FLAG) Positive Control Lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

**Clonality**

Polyclonal

**Dilution Range**

WB (1:1000)

**Expiration Date**

12 months from date of receipt.