



Product Datasheet

CD122 antibody (orb414539)


Description	Mouse monoclonal antibody to CD122
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP
Immunogen	TL-Mor cell line
Target	CD122
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P14784
Entrez	3560
Dilution Range	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD122
positive CD56
...



Flow
cytometry
surface
staining
pattern ...