

Product Datasheet

NEDD4L Antibody, FITC conjugated (orb415215)

Catalog Number	orb415215
Category	Antibodies
Description	This NEDD4L Antibody, FITC conjugated is an FITC polyclonal product. It tagets NEDD4L using a Recombinant Human E3 ubiquitin-protein ligase NEDD4-like protein (289-494AA) as the immunogen. Purification: >95%, Protein G purified.
Target	NEDD4L
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human E3 ubiquitin-protein ligase NEDD4-like protein (289-494AA)
UniProt ID	Q96PU5
Dilution range	ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558