



Product Datasheet

IFNGR2 antibody (orb415521)


| | |
|----------------------------|--|
| Description | IFNGR2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC |
| Immunogen | Recombinant Human Interferon gamma receptor 2 protein (28-205AA) |
| Target | IFNGR2 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P38484 |
| Dilution Range | IHC-P:1:20-200, IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemistry of paraffin-embedde...



Immunofluorescence staining of PC-3 cell...



Immunohistochemistry of paraffin-embedde...