

Product Datasheet

NDUFS2 antibody (orb415641)

Description	NDUFS2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] iron-sulfur protein 2, mitochondrial protein (297-427AA)
Target	NDUFS2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75306
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.