

Product Datasheet

MKRN1 antibody (HRP) (orb416170)

Description	MKRN1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human E3 ubiquitin-protein ligase makorin-1 protein (109-209AA)
Target	MKRN1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UHC7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.