

Product Datasheet

ARHGAP10 antibody (HRP) (orb41650)

Description	ARHGAP10 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Rho GTPase-activating protein 10 protein (542-744AA)
Target	ARHGAP10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A1A4S6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.