



Product Datasheet

HIST1H4A (Ab-8) antibody (orb416507)


Description	HIST1H4A (Ab-8) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, IP, WB
Immunogen	Peptide sequence around site of Lys (8) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot analysis of Rat liver tiss...



Immunoprecipitation analysis of NIH/3T3

...