



Product Datasheet

Tri-methyl-HIST1H3A (K4) antibody (orb416548)


Description	Tri-methyl-HIST1H3A (K4) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Peptide sequence around site of Tri-methyl-Lys (4) derived from Human Histone H3.1
Target	HIST1H3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431
Dilution Range	WB:1:500-5000, ICC:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cells...



Western blot analysis of Jurkat cell lysates...



Immunofluorescent analysis of HeLa cells...