



Product Datasheet

HIST1H1E (phospho-S35) antibody (orb416595)

Description	HIST1H1E (phospho-S35) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Peptide sequence around site of Phospho-Ser (35) derived from Human Histone H1.4
Target	HIST1H1E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10412
Dilution Range	ICC:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cel...



Immunofluorescent analysis of HeLa cells...