



Product Datasheet

H2AFZ (Ab-4) antibody (orb416648)


Description	H2AFZ (Ab-4) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, WB
Immunogen	Peptide sequence around site of Lys (4) derived from Human Histone H2A.Z
Target	H2AFZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C055
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot analysis of HepG2 whole ce...



Immunofluorescence staining of Hela cell...