



Product Datasheet

BUB1B antibody (orb41677)

Description	BUB1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 beta protein (230-394AA)
Target	BUB1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60566
Dilution Range	IHC-P:1:20-200, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human en...



Immunofluorescence analysis of HepG2 cel...