

Product Datasheet

Beta galactosidase antibody (HRP) (orb416772)

Description	Beta galactosidase antibody (HRP)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Beta-galactosidase protein (2-1024AA)
Target	lacZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00722
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.