

Product Datasheet

TNIP1 antibody (Biotin) (orb417098)

Description	TNIP1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human TNFAIP3-interacting protein 1 protein (526-636AA)
Target	TNIP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15025
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.