



Product Datasheet

PAFAH1B3 antibody (orb417099)


Description	PAFAH1B3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Platelet-activating factor acetylhydrolase IB subunit gamma protein (103-187AA)
Target	PAFAH1B3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15102
Dilution Range	IHC-P:1:500-1000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of A549 cell...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...